



People.Health.Care.

Velpeau® Manus Comfort **NEW**

Active wrist support.

COMFORT



Velpeau Manus Comfort stabilises the wrist in a functionally neutral position, provides compression and reduces swelling.

- provides compression through the combination of the knitted fabric, the crescent-shaped silicone insert and the stabilising strap
- knitted three-dimensional anatomical shape
- special seamless knitted fabric in the sensitive thumb area





Individually adjustable

- the splint can be thermoplastically moulded as required



Physiological fit

- three-dimensional shaped knitted support
- recesses in the strap and insert prevent excessive pressure on bony protuberances



High wearing comfort

- designed to be comfortable due to the microfibre materials
- specially knitted fabric in the thumb area to selectively avoid pressure points

Velpeau® Manus Comfort

Active wrist support

Wrist circumference



Size	Wrist circumference (cm)	REF	Shipping units (pcs.)
Left hand version			
T0 G	13.0-14.5	34 518	1 / 64
T1 G	14.5-16.5	34 519	1 / 64
T2 G	16.5-19.0	34 520	1 / 64
T3 G	19.0-23.0	34 521	1 / 64
Right hand version			
T0 G	13.0-14.5	34 484	1 / 64
T1 G	14.5-16.5	34 485	1 / 64
T2 G	16.5-19.0	34 486	1 / 64
T3 G	19.0-23.0	34 487	1 / 64



Indications

- activated osteoarthritis or arthritis of the wrist
- moderately severe to severe distortions (sprains) in post-trauma patients
- tendinitis in the hand region (tendovaginitis)
- instability and functional overloading in the wrist area
- follow-up treatment for wrist fractures, after surgery or after wearing a plaster cast
- stage I complex regional pain syndrome (CRPS)

Colour

taupe, white