

Velpeau® Freecast® Adjust Classic **NEW**

Orthosis for stabilisation of the ankle joint.

CLASSIC



Velpeau Freecast Adjust Classic stabilises the ankle in one plane based on the principle of the standard U splint.

- anatomically shaped shells
- individually adjustable
- one size fits all, for use on right or left ankle





Safe to use

- easy to fasten with hook-and-loop fasteners
- visual application instructions on the product



One size fits all

- individually adjustable to the patient's anatomy
- easy to stock and store



High wearing comfort

- soft, flexible area for the malleoli
- internal padding with a mesh cover for good ventilation

Velpeau® Freecast® Adjust Classic

Stabilising ankle support

Size	Shoe size	REF	Shipping units (pcs.)
T1	35-45	34 007	1 / 30



Indications

- ligament injuries to the upper ankle joint (grade 1 – 2)
- activated osteoarthritis, arthritis
- chronic ligament instability
- perimalleolar swelling, joint effusion, intra-articular irritation

Colour

black, green